


# LIBERTY

 **EPS 2.0**  
makes the  
difference



**MEDICENTUR**

technology upgrade

# Technology Comparison




Enhanced and more homogenous intermediate light intensity results in better visual quality


Effective light intensity pattern among the three focal points on PSF (Point Spread Function) Maps in ISO Model Eye.

The lightness on the optical axis correlates with image sharpness and contrast sensitivity at the given distance. Light detected away from the optical axis could be observed as dysphotopsias.


**Liberty  
EPS 2.0**




**Liberty  
EPS 1.0**



**PanOptix**



**FineVision**



Log of  
relative  
intensity[-]




# Technology Comparison



Less out of focus light results in less diffractive dysphotopsias

Effective light intensity pattern at the far focal point on front PSF Maps in ISO Model Eye simulated in scotopic conditions.

The lightness on the center correlates with image sharpness and contrast sensitivity at the given distance. Light detected away from the center could be observed as dysphotopsias.



# Technology Comparison




## Axial MTF in ISO Model Eye: improved intermediate in EPS 2.0


### Axial MTF map

Axial Modulation Transfer Function (MTF) maps show the combination of 2D MTF and Through Focus MTF characteristics at a given aperture, providing a **visual representation of the contrast sensitivity distribution** and vision quality at different distances.

The images show MTF intensities (colouring) through different spatial frequencies (y) regarding the defocus power (x). MTF distribution is measured on a 3 mm aperture, the average pupil diameter in photopic conditions.

Compared to the EPS1.0 technology, **the EPS 2.0 provides a remarkably higher intensity in the intermediate region**, especially at 15 cpd spatial frequency. In comparison with the two main market leaders (PanOptix and FineVision), **EPS 2.0 shows a more homogenous intermediate and more intense near MTF, resulting in superior visual quality at all distances.**

The [QR code directs you](#)  to a simulation showing the MTF distribution change at ascending aperture widths. The MTF intensity dominance changes with the aperture between distances, as pupil diameter changes in different environments. Near power is more intense in photopic conditions (for reading), while far vision is most dominant in scotopic conditions, enabling the possibility of night driving.



The light intensity measurements support the MTF simulation: as the aperture changes, the light distribution among the focal points changes accordingly.





# Technology Comparison



Extended effective MTF region with higher MTF values in the intermediate range results in higher contrast sensitivity and better visual acuity

Axial MTF in ISO Model Eye: improved intermediate in EPS 2.0.



# Clinical Outcomes




Confirm improved visual performance

Binocular visual acuity and contrast sensitivity defocus curves confirm **excellent visual acuity and visual quality** throughout the entire defocus range.<sup>1</sup>




EPS 2.0 Intermediate vision improved significantly compared to the previous technology without compromising distance and near vision.<sup>1</sup>



a

Improved  
Intermediate  
Vision  
in the range  
from -1 to -2 D



b

Provides Excellent  
Visual Performance  
without any loss in  
far and near vision

Both defocus curves  
were obtained 3 months  
postoperatively under  
photopic light conditions  
(85 cd/m<sup>2</sup>).

1. Fernández J. EPS – Pushing the limits in shifting the wavefront.  
Presented during the Faco Elche congress in 2022, Elche, Alicante, Spain.

Measurements were performed using the Qvision Multifocal Lens Analyzer. (Fernández, J. et al. Fast Measure of Visual Acuity and Contrast Sensitivity Defocus Curves with an iPad Application. Open Ophthalmol J 13, 15–22 (2019).

# Clinical Measurements:



## Low LDI (Light Distortion Index) values: high quality vision



The subjective perception of light distortion is an indicator of visual quality. The lower the light distortion index (LDI), the better the visual quality is.



The Liberty IOLs have lower LDI values compared to other multifocal IOLs on the market. The improved EPS 2.0 models' LDI is even lower than a monofocal lens.

Author	Year	IOL	LDI	
Short term			Monocular (%)	Binocular (%)
Brito P et al. <sup>2</sup>	2015	Tecnis ZCB00 ( <b>Monofocal</b> )	24	15
		AT Lisa 839M (Tri)	47	29
		AT Lisa 909MP (Bi-Tor)	54	40
Alió et al. <sup>3</sup>	2018	AcrySof IQ Panoptix	37	24
Oliveira RF et al. <sup>4</sup>	2018	FineVision	33	23
Escandón García et al. <sup>5</sup>	2018	FineVision		29
		AcrySof IQ Panoptix		26
		Tecnis Symphony		35
Vargas V et al. <sup>6</sup>	2020	MPlus +1.5 / +3.0 D	21 / 27	
Fernández et al. (Publication pending)	3 month 12 month	<b>Liberty EPS 2.0 (640PM)</b>	<b>14 ± 7 15 ± 8</b>	<b>13 ± 8</b>

The light-transmission analyser is a useful tool in the postoperative evaluation of the quality of vision conferred by multifocal IOLs. The software calculates the distortion index (DI) as the ratio of the area of points missed by the subject and the total area explored and is expressed as a percentage. The deviation of the obtained polygonal shape from the best-fit circle fit is called the best-fit circle irregularity (BFCI).<sup>1-2</sup>



### Size & Shape Parameters

- Light Distortion/Disturbance Index - LDI [%]
- Best Fit Circle Radius BFC\_Radius [%]

### Regularity Parameters

- LD Irregularity BFC\_Irregularity [mm]
- SD LD Irregularity BFC\_IrregSD [mm]

2. Brito P et al. J Cataract Refract Surg. 2015;41(3):613-22
3. Alió J et al. Eur J Ophthalmol. 2018;28(4):419-24
4. Oliveira RF et al. EJO. 2018;1-8
5. Escandón-García S et al. J Ophthalmol. 2018;2018
6. Vargas V et al. J Refract Surg. 2020;36(12):796-803



## **Medicontur Medical Engineering Ltd.**

### **Headquarters:**

Herceghalmi Road 1  
2072 Zsámbék  
Budapest / Hungary

### **International Sales Office:**

Chemin du Champ-des-Filles 36.  
1228 Plan-les-Ouates  
Switzerland / Geneva

### **More information:**

[www.medicontur.com](http://www.medicontur.com)

