



ASICO®

IOL Cutter & Scissors




AE-4296, ASICO IOL Cutter

Designed to securely cut and remove silicone and/or acrylic IOLs. Atraumatic shape on the front end of Features a cavity for holding the lens. the tip for easy insertion through a 3mm wound.



AE-5437, Shimowake IOL Scissors

Designed to cut an acrylic/silicone lens. Tips also feature a rounded ball for added safety, designed to prevent any potential tears on the capsule or iris. Angled tips with serrated edge.



AE-5438, Akahoshi IOL Scissors

Designed tips to cut an acrylic/silicone lens and the nucleus. For added safety, the tips of the scissors have a rounded ball to prevent any potential tears on the capsule and/or the iris.

The more you know...
ASICO®
Today's Precision... Tomorrow's Vision™

ASICO LLC
26 Plaza Drive • Westmont, IL 60559 • USA
Phone: 630-986-8032 • Fax: 630-986-0065
Customer Service: 1-800-628-2879
www.asico.com • E-mail: info@asico.com
©2012 ASICO LLC

