

Technical Data Sheet – IVT Procedure Pack

Product code: GBSPUI06-05

Case quantity: 50

Content:

- 1 Wrap crepe 90 x 90cm
- 1 Tray 3 compartment 13 x 20cm
- 5 Gauze swab 10 x 10cm 12ply
- 2 Cotton bud small synthetic grip 7.5cm
- 1 Speculum Barraquer temporal solid blade
- 1 Double ended marker plastic 3.5/4mm
- 1 Transparent drape 50 x 40cm aperture 7x10cm and perforated incision film

Storage: Store in a dry and cool place, away from intense sources of heat.
Shelf life: 4 years, from date of manufacture.
Delivery time: 1-2 days after placing order.
Minimum Order: 1 case.

Manufacturing: Products are manufactured in Germany. The quality system of manufacturing sites is ISO 13485 compliant

Regulatory info: All our products have a CE mark and comply with the requirements of the Medical Device Directive 93/42/EEC Class: IIa



Increasing volume of IVTs? Injection rooms a limited resource?
Get in touch for more info on Operio Toul Mobile and IVT pack offers to utilise additional clinic rooms and improve efficiency!