

Technical Data Sheet

GBSPUI06-05 – IVT Procedure Pack

Content:	1	Wrap crepe 90 x 90cm
	1	Tray 3 compartment 13 x 20cm
	5	Gauze swab 10 x 10cm 12ply
	2	Cotton bud small synthetic grip 7.5cm
	1	Speculum Barraquer temporal solid blade
	1	Double ended marker plastic 3.5/4mm
	1	Transparent patient drape 50 x 40cm aperture 7x10cm and perforated incision film



Main materials: Wrap: semicreped paper
Drape: PE-foil, latex-free
Does not contain DEHP



Manufacturing: Products are manufactured in Germany.
The quality system of manufacturing sites is ISO 13485 compliant

Regulatory info: All our products have a CE mark and comply with the requirements of the Medical Device Directive 93/42/EEC Class: IIa

Storage: Store in a dry and cool place, away from intense sources of heat.
Shelf life: 4 years, from date of manufacture.
Delivery time: 2 days after placing order.
Minimum Order: 1 case.



Increasing volume of IVTs? Injection rooms a limited resource?
Get in touch for more info on Operio Toul Mobile and IVT pack offers to utilise additional clinic rooms and improve efficiency!

