


SWISS COAXIAL SOFT TIP I/A



SWISS COAXIAL SOFT TIP I/A



The gentle alternative

The SOFT TIP I/A, made of soft silicone, allows gentle contact with the sensitive capsular bag tissue. The textured handpiece grip surface offers comfortable handling and user control.

SWISS COAXIAL SOFT TIP I/A system, with silicone tip (single-use / sterile) consisting of: Coaxial I/A instrument (box of 10)

Item Code	Size	Incision Size*	Aspiration Cannula	Aspiration Port	Aspiration Head
SMS390S	2IG / 0.8mm	2.2mm	angled 45°	0.4mm	silicone tip

* Incision sizes are indications only based on average experiences. The individually required incision sizes may vary depending on surgical techniques



instruction for use



Medicel AG, Dornierstrasse 11, 9423 Altenrhein, +41 71 727 10 50, medicel.com



digital brochure

made in switzerland